

HCV NS5 Genotype-1a Recombinant Protein



RPPB5467

Product Information Protein Information

Product SKU:

RPPB5467

Protein description:

The E.coli derived 38 kDa recombinant protein fused to a His tag contains the HCV NS5 Genotype1a immunodominant regions, amino acids 2212-2313.

Formulation:

1.5M urea, 25 mM Tris-HCl, pH-8, 0.2% Triton-X & 50% Glycerol.

Purity:

HCV NS5 Genotype-1a protein is >95% pure as determined by 10% PAGE (coomassie staining).

Stability:

HCV NS5 Genotype-1a although stable at 4°C for 1 week, should be stored below -18°C. Please prevent freeze thaw cycles.