

# HCV NS5 Genotype-1b Recombinant Protein



RPPB5468

---

## Product Information    Protein Information

**Product SKU:**

RPPB5468

**Host:**

Escherichia Coli.

**Protein description:**

The E.coli derived recombinant protein contains the HCV NS5 Genotype1b immunodominant regions, amino acids 2212-2313. The protein is fused to a GST tag at N-terminus.

**Appearance:**

Sterile Filtered clear solution.

**Formulation:**

1.5M urea, 25mM Tris-HCl pH 8.0, 0.2% Triton-X and 50% Glycerol.

**Purity:**

HCV NS5 Genotype-1b protein is >95% pure as determined by SDS- PAGE (coomassie staining).

**Stability:**

HCV NS5 Genotype-1b although stable at 4°C for 1 week, should be stored below -18°C. Please prevent freeze thaw cycles.