

# HCV NS5 Genotype-2 Recombinant Protein



RPPB5469

---

## Product Information    Protein Information

**Product SKU:**

RPPB5469

**Protein description:**

The E.coli derived recombinant protein contains the HCV NS5 Genotype2 immunodominant regions.

**Formulation:**

25mM Tris and 1mM EDTA, 1.5M urea & 50% glycerol.

**Purity:**

HCV NS5 Genotype-2 protein is >95% pure as determined by 10% PAGE (coomassie staining).

**Stability:**

HCV NS5 Genotype-2 although stable at 4°C for 1 week, should be stored below -18°C. Please prevent freeze thaw cycles.