

HCV NS5 Genotype-2b Recombinant Protein



RPPB5471

Product Information Protein Information

Product SKU:

RPPB5471

Host:

Escherichia Coli.

Protein description:

The E.coli derived recombinant protein contains the HCV NS5 Genotype 2b immunodominant regions, amino acids 2212-2313. The protein is fused to a GST tag at N-Terminus.

Formulation:

1.5M urea, 25 mM Tris-HCl, pH-8, 0.2% Triton-X & 50% Glycerol.

Purity:

HCV NS5 Genotype-2b protein is >95% pure as determined by 10% PAGE (coomassie staining).

Stability:

HCV NS5 Genotype-2b although stable at 4°C for 1 week, should be stored below -18°C. Please prevent freeze thaw cycles.