

HCV NS5 Genotype-3 Recombinant Protein



RPPB5472

Product Information Protein Information

Product SKU:

RPPB5472

Protein description:

The E.coli derived recombinant protein fused to a GST tag contains the HCV NS5 Genotype 3 immunodominant regions.

Formulation:

25mM Tris-HCl, 1mM EDTA, 1.5M urea and 50% glycerol.

Purity:

HCV NS5 Genotype-3 protein is >95% pure as determined by 10% PAGE (coomassie staining).

Stability:

HCV NS5 Genotype-3 although stable at 4°C for 1 week, should be stored below -18°C. Please prevent freeze thaw cycles.