

HCV NS5 Genotype-4 Recombinant Protein



RPPB5475

Product Information Protein Information

Product SKU:

RPPB5475

Protein description:

The E.coli derived recombinant protein contains the HCV NS5 Genotype 4 immunodominant regions, amino acids 2212-2313. The protein is fused with GST tag at N-terminus.

Formulation:

1.5M urea, 25mM Tris-HCl, pH-8, 0.2% Triton-X & 50% Glycerol.

Purity:

HCV NS5 Genotype-4 protein is >95% pure as determined by 10% PAGE (coomassie staining).

Stability:

HCV NS5 Genotype-4 although stable at 4°C for 1 week, should be stored below -18°C. Please prevent freeze thaw cycles.