

HCV NS5 Genotype-6 Recombinant Protein



RPPB5477

Product Information Protein Information

Product SKU:

RPPB5477

Protein description:

The E.coli derived recombinant protein contains the HCV NS5 Genotype 6 immunodominant regions.

Formulation:

25mM Tris-HCl, 1mM EDTA, 1.5M urea and 50% glycerol.

Purity:

HCV NS5 Genotype-6 protein is >95% pure as determined by 10% PAGE (coomassie staining).

Stability:

HCV NS5 Genotype-6 although stable at 4°C for 1 week, should be stored below -18°C. Please prevent freeze thaw cycles.