

# HCV NS5 Recombinant Protein



RPPB5479

---

## Product Information    Protein Information

**Product SKU:**

RPPB5479

**Protein description:**

The E.coli derived recombinant protein contains the HCV NS5a immunodominant regions, amino acids 2061-2302. The protein is fused with GST at N-terminus.

**Formulation:**

1.5M Urea, 25mM Tris-HCl pH 8, 50% glycerol and 0.2% Triton-X.

**Purity:**

HCV-NS5 protein is >95% pure as determined by 10% PAGE (coomassie staining).

**Stability:**

HCV NS5 although stable at 4°C for 1 week, should be stored below -18°C. Please prevent freeze thaw cycles.