

HCV NS5B Recombinant Protein



RPPB5481

Product Information Protein Information

Product SKU:

RPPB5481

Protein description:

The E.coli derived recombinant protein contains the full length of the NS5B immunodominant regions Genotype 1B and fused to a 6xHis tag at C-Terminus.

Formulation:

1.5M urea, 25mM Tris-HCl pH 8.0, 0.2% Triton-X and 50% Glycerol.

Purity:

Protein is >95% pure as determined by SDS-PAGE.

Stability:

HCV NS5B although stable at 4°C for 1 week, should be stored below -18°C. Please prevent freeze thaw cycles.