

HCV NS5 Recombinant Protein



RPPB5483

Product Information Protein Information

Product SKU:

RPPB5483

Protein description:

The E.coli derived Biotin Labeled recombinant protein contains the HCV NS5 immunodominant regions. HCV NS5 antigen (recombinant) a.a 2061 to a.a 2302 of HCV polyprotein. The protein is fused to a GST tag at N-terminus.

Formulation:

1.5M urea, 25mM Tris-HCl pH 8.0, 0.2% Triton-X and 50% Glycerol.

Purity:

HCV NS5 Biotin protein is >95% pure as determined by 10% PAGE (coomassie staining).

Stability:

HCV NS5 Biotin although stable at 4°C for 1 week, should be stored below -18°C. Please prevent freeze thaw cycles.