

HCV NS5 Recombinant Protein



RPPB5484

Product Information Protein Information

Product SKU:

RPPB5484

Protein description:

The E.coli derived Biotin Labeled recombinant protein contains the HCV NS5 immunodominant regions.

Formulation:

50mM Tris pH-8 & 5mM EDTA.

Purity:

HCV-NS5 HRP protein is >95% pure as determined by 10% PAGE (coomassie staining).

Stability:

HCV NS5 HRP although stable at 4°C for 1 week, should be stored below -18°C. Please prevent freeze thaw cycles.