

HCV 4th Generation Recombinant Protein



RPPB5485

Product Information Protein Information

Product SKU:

RPPB5485

Host:

Escherichia Coli.

Protein description:

The E.Coli derived HCV fourth generation antigen recombinant protein contains medium size core (55aa), NS3 (226aa), NS4 (3 epitopes) & NS5 (3 epitopes).HCV 4th Generation migrates at 46kDa.

Formulation:

Phosphate saline buffer with 25mM arginine.

Purity:

Protein is >95% pure as determined by 12% PAGE (coomassie staining).

Stability:

HCV 4th Generation although stable at 4°C for 1 week, should be stored below -18°C. Please prevent freeze thaw cycles.