

HAV VP1-P2A (722-830 a.a.) Recombinant Protein



RPPB5490

Product Information Protein Information

Product SKU:

RPPB5490

Protein description:

The E.Coli derived 51.2 kDa recombinant protein fused to a GST tag contains the VP1-P2A immunodominant regions, amino acids 722-830.

Formulation:

10mM CBB, pH9.6, 0.1% SDS and 50% glycerol.

Purity:

HAV VP1-P2A protein is >90% pure as determined by 10% PAGE (coomassie staining).

Stability:

HAV VP1-P2A protein although stable at 4°C for 1 week, should be stored below -18°C. Please prevent freeze thaw cycles.