

# HAV VP1-P2A (669-782 a.a) Recombinant Protein



RPPB5491

---

## Product Information    Protein Information

**Product SKU:**

RPPB5491

**Host:**

Escherichia Coli.

**Protein description:**

The E.Coli derived recombinant protein contains the VP1-P2A immunodominant regions, amino acids 669-782.

**Formulation:**

10mM Tris-HCl, pH 9.6, 1.5M urea, 50% glycerol.

**Purity:**

HAV VP1-P2A protein is >90% pure as determined by 10% PAGE (coomassie staining).

**Stability:**

HAV VP1-P2A protein although stable at 4°C for 1 week, should be stored below -18°C. Please prevent freeze thaw cycles.