

# HAV P2C Recombinant Protein



RPPB5492

---

## Product Information    Protein Information

**Product SKU:**

RPPB5492

**Protein description:**

The E.Coli derived recombinant 38.5 kDa protein contains the P2C immunodominant regions, amino acids 1121-1234.

**Formulation:**

10mM CBB, pH9.6, 0.1% SDS and 50% glycerol.

**Purity:**

Protein is >90% pure as determined by 10% PAGE (coomassie staining).

**Stability:**

HAV P2C protein although stable at 4°C for 1 week, should be stored below -18°C. Please prevent freeze thaw cycles.