

HAV P2C-P3B Recombinant Protein



RPPB5493

Product Information Protein Information

Product SKU:

RPPB5493

Host:

Escherichia Coli.

Protein description:

The E.Coli derived 53.7 kDa recombinant protein fused to a GST tag contains the P2C-P3B immunodominant regions, amino acids 1492-1606.

Formulation:

10mM Tris-HCl, pH 9.6, 0.1% SDS and 50% glycerol.

Purity:

HAV P2C-P3B protein is >90% pure as determined by 10% PAGE (coomassie staining).

Stability:

HAV P2C-P3B Protein although stable at 4°C for 1 week, should be stored below -18°C. Please prevent freeze thaw cycles