

HAV P2C-P3A Recombinant Protein



RPPB5494

Product Information Protein Information

Product SKU:

RPPB5494

Host:

Escherichia Coli.

Protein description:

The E.Coli derived recombinant protein contains the P2C-P3A immunodominant regions, amino acids 1392-1521.

Formulation:

10mM CBB, pH9.6, 0.1% SDS and 50% glycerol.

Purity:

HAV P2C-P3A protein is >90% pure as determined by 10% PAGE (coomassie staining).

Stability:

HAV P2C-P3A protein although stable at 4°C for 1 week, should be stored below -18°C. Please prevent freeze thaw cycles.