

HAV P3C Recombinant Protein



RPPB5495

Product Information Protein Information

Product SKU:

RPPB5495

Host:

Escherichia Coli.

Protein description:

The E.Coli derived 40 kDa recombinant protein contains the P3C immunodominant regions, amino acids 1643-1743.

Formulation:

10mM CBB, pH9.6, 0.1% SDS and 50% glycerol.

Purity:

HAV P3C protein is >90% pure as determined by 10% PAGE(coomassie staining).

Stability:

HAV P3C protein although stable at 4°C for 1 week, should be stored below -18°C. Please prevent freeze thaw cycles.