

# HAV VP3 Recombinant Protein



RPPB5496

---

## Product Information    Protein Information

**Product SKU:**

RPPB5496

**Host:**

Escherichia Coli.

**Protein description:**

The E.coli derived 38 kDa recombinant protein contains the VP3 immunodominant regions, amino acids 304-415.

**Formulation:**

10mM CBB, pH9.6, 0.1% SDS and 50% glycerol.

**Purity:**

HAV VP3 protein is >90% pure as determined by 10% PAGE (coomassie staining).

**Stability:**

HAV VP3 Protein although stable at 4°C for 1 week, should be stored below -18°C. Please prevent freeze thaw cycles.