

HAV VP4-VP2 Recombinant Protein



RPPB5497

Product Information Protein Information

Product SKU:

RPPB5497

Host:

Escherichia Coli.

Protein description:

The E.Coli derived 44 kDa recombinant protein fused to a GST tag contains the VP4-VP2 immunodominant regions, amino acids 55-164.

Formulation:

10mM CBB, pH9.6, 0.1% SDS and 50% glycerol.

Purity:

HAV VP4-VP2 protein is >90% pure as determined by 10% PAGE (coomassie staining).

Stability:

HAV VP4-VP2 protein although stable at 4°C for 1 week, should be stored below -18°C. Please prevent freeze thaw cycles.