

HBV core (1-186 a.a.) Recombinant Protein



RPPB5503

Product Information Protein Information

Product SKU:

RPPB5503

Protein description:

The E.coli derived recombinant protein contains the HBV core immunodominant region amino acids 1-186, and fused to a His tag at N- terminus.

Formulation:

25mM Tris-HCl pH-8.0, 1.5mM EDTA, 1.5mM Urea and 50% glycerol.

Purity:

HBV Core protein is >90% pure as determined by 10% PAGE (Coomassie staining).

Stability:

HBV Core protein although stable at 4°C for 1 week, should be stored below -18°C. Please prevent freeze thaw cycles.