

HBV core (1-183 a.a.) Recombinant Protein



RPPB5504

Product Information Protein Information

Product SKU:

RPPB5504

Protein description:

The E.Coli derived 18 kDa recombinant protein contains the HBV core ayw immunodominant region, amino acids 1-183.

Formulation:

1x PBS (pH 7.2), 75 mM NaCl and 20% glycerol.

Purity:

HBV Core protein is >95% pure as determined by SDS-PAGE.

Stability:

HBV Core protein although stable at 4°C for 1 week, should be stored below -18°C. Please prevent freeze thaw cycles.

Amino Acid Sequence:

MDIDPYKEFG ASVELLSFLP SDFFPSIRD LDTASALYRE ALESPEHCSPHHTALRQAIL CWGELMNLAT
WVGSNLEDPA SRELVSYSVN VNMGLKFRQL LWFHVSCLTF GRETVLEYLV SFGVWIRTPP AYRPPNAPIL
STLPETTIVR RRGRSPRRRT PSPRRRRSQS PRRRRSQSRE SQ.