

HBV core delta Recombinant Protein



RPPB5505

Product Information Protein Information

Product SKU:

RPPB5505

Protein description:

The E.coli derived 14kDa recombinant protein contains the HBV core delta ayw immunodominant region amino acids 1-144.

Formulation:

7.5mM Phosphate buffer pH-7.2, 75mM NaCl, & 50% glycerol.

Purity:

HBV Core Delta protein is >95% pure as determined by SDS-PAGE (coomassie staining).

Stability:

HBV Core Delta protein although stable at 4°C for 1 week, should be stored below -18°C. Please prevent freeze thaw cycles.

Amino Acid Sequence:

MDIDPYKEFG ATVELLSFLP SDFFPSVRDL LDNASALYRE ALESPEHCSP HHTALRQAIL CWGELMTLAT
WVGGNLEDPI SRDLVSVYVN TNMGLKFRQL LWFHISCLTF GRETVIEYLV SFGVWIRTPP AYRPPNAPIL S