

HBV HBe Recombinant Protein



RPPB5510

Product Information Protein Information

Product SKU:

RPPB5510

Protein description:

The E.Coli derived recombinant protein contains the HBV HBe adw immunodominant region is fused to a GST tag and the Mw is 43,897.6 Dalton.

Formulation:

(1mg/ml) 50mM Tris-HCl pH 8.5, 5mM DTT, 10mM BME and 8M Urea.

Purity:

HBV HBe protein is >95% pure as determined by 10% PAGE (coomassie staining).

Stability:

HBV HBe protein although stable at 4°C for 1 week, should be stored below -18°C. Please prevent freeze thaw cycles.