

HBV Pre-S Recombinant Protein



RPPB5511

Product Information Protein Information

Product SKU:

RPPB5511

Host:

Escherichia Coli.

Protein description:

Recombinant Hepatitis B Virus Pre-S produced in E.coli is a full length Pre-S of HBV subtype adw2, containing 174 amino acids, having a molecular weight of about 25kDa including a vector sequence and linked with a 6xHis tag at C- terminus.

Appearance:

Sterile Filtered solution.

Formulation:

Sterile Filtered solution containing PBS, 1M urea and 0.05% sodium azide.

Purity:

Protein is >95% pure as determined by 12% PAGE (coomassie staining).

Stability:

HBV Pre-S although stable at 4°C for 1 week, should be stored below -18°C. Please prevent freeze thaw cycles.