

HTNV Recombinant Protein



RPPB5512

Product Information Protein Information

Product SKU:

RPPB5512

Host:

E.Coli

Protein description:

The Hantavirus nucleocapsid (N) fusion protein is expressed in E. coli and is fused to 6xHis tag. Hantavirus nucleocapsid (N) is conserved among different strains, and is used to test the specific IgM and IgG to Hantavirus.

Appearance:

Sterile filtered colorless solution.

Synonyms:

HTNV, Hantaan Virus, Hantanvirus.

Formulation:

HTNV is formulated in 1xPBS pH-7.4 and 0.05 Sodium azide.

Purity:

HTNV protein is >95% pure as determined by 12% PAGE (coomassie staining).

Stability:

Upon arrival, Store at -20°C. Please avoid multiple freeze/thaw cycles.

Amino Acid Sequence:

GSMATMEELQREINAHEGQLVIARQKVRDAEKQYEKDPDELNKRTLTDRE
GVAVSIQAKIDELKRQLADRIATGKNLGKEQDPTGVEPGDHLKERSMLSY
GNVLDLNHLDIDEPTGQTADWLSIVVYVD