

FeLV Recombinant Protein



RPPB5513

Product Information Protein Information

Product SKU:

RPPB5513

Host:

Escherichia Coli.

Protein description:

The E.Coli truncated derived recombinant protein contains the p27 immunodominant regions, containing 248 amino acids (271-519) and having a Mw of 31kda as visualized by SDS-PAGE. Recombinant Feline Leukemia Virus p27 is fused to a His Tag and purified by conventional chromatography techniques.

Appearance:

Sterile Filtered clear solution.

Formulation:

1x PBS

Purity:

Protein is >95% pure as determined by 10% PAGE (coomassie staining).

Stability:

FeLV p27 although stable at 4°C for 1 week, should be stored below -18°C. Please prevent freeze thaw cycles.