

# HEV Mosaic Recombinant Protein



RPPB5515

---

## Product Information   Protein Information

**Product SKU:**

RPPB5515

**Protein description:**

The e.coli derived HEV Mosaic protein contains 12 HEV immunodominant regions from ORF2 and ORF3 having a molecular mass of 38.5 kDa.

**Formulation:**

25mM Tris pH 8.0, 1mM EDTA, 0.5M urea and 50% glycerol.

**Purity:**

HEV Mosaic protein is >95% pure as determined by 10% PAGE (coomassie staining).

**Stability:**

HEV Mosaic protein although stable at 4°C for 1 week, should be stored below -18°C. Please prevent freeze thaw cycles.