

HEV ORF2 (633-659 a.a.) Recombinant Protein



RPPB5516

Product Information Protein Information

Product SKU:

RPPB5516

Protein description:

The E.coli derived HEV protein genotype-2 is fused with beta-galactosidase at the N-Terminus and contains the HEV immunodominant ORF2 633-659 a.a.

Formulation:

20mM Tris-HCl pH-8, 8M urea and 10mM B-ME.

Purity:

HEV ORF2 Protein is >95% pure as determined by 10% PAGE (coomassie staining).

Stability:

HEV ORF2 protein although stable at 4°C for 1 week, should be stored below -18°C. Please prevent freeze thaw cycles.