

HEV ORF2 (403-461 a.a.) Recombinant Protein



RPPB5517

Product Information Protein Information

Product SKU:

RPPB5517

Protein description:

The e.coli derived HEV protein is fused with beta-galactosidase at the N-Terminus and contains the HEV immunodominant regions from ORF2 403-461 a.a.

Formulation:

20mM Tris-HCl pH-8, 8M urea and 10mM B-ME.

Purity:

HEV ORF2 Protein is >95% pure as determined by 10% PAGE (coomassie staining).

Stability:

HEV ORF2 protein although stable at 4°C for 1 week, should be stored below -18°C. Please prevent freeze thaw cycles.