

HEV ORF2 (452-617 a.a.) Recombinant Protein



RPPB5518

Product Information Protein Information

Product SKU:

RPPB5518

Protein description:

The E.Coli derived HEV protein contains the HEV immunodominant regions from ORF2 452-617 a.a. having a Mw. of 44.3 kDa. HEV ORF2 protein is fused to a 26kDa GST tag and purified by standard chromatography techniques.

Formulation:

25mM Tris-HCl, 1mM EDTA, 0.5M urea and 50% glycerol.

Purity:

HEV ORF2 protein is >95% pure as determined by 10% PAGE (coomassie staining).

Stability:

HEV ORF2 protein although stable at 4°C for 1 week, should be stored below -18°C. Please prevent freeze thaw cycles.