

HEV ORF3 Recombinant Protein



RPPB5519

Product Information Protein Information

Product SKU:

RPPB5519

Protein description:

The e.coli derived HEV protein is fused with beta-galactosidase at the N-Terminus and contains the HEV immunodominant ORF3 92-123 a.a.

Formulation:

20mM Tris-HCl, pH-8, 10mM B-ME and 8M urea.

Purity:

Greater than 95.0% as determined by Analysis by SDS-PAGE.

Stability:

HEV ORF3 protein although stable at 4°C for 1 week, should be stored below -18°C. Please prevent freeze thaw cycles.