

HBsAg adr CHO Recombinant Protein



RPPB5522

Product Information Protein Information

Product SKU:

RPPB5522

Protein description:

Recombinant HBsAg adr produced in CHO cells, the molecular weight is approximately 23-27kDa.

Accession:

P30019

Appearance:

Sterile Filtered clear solution.

Host:

CHO.

Formulation:

Sterile Filtered solution containing 20mM PB and 154mM NaCl.

Purity:

Greater than 95.0% as determined by RP-HPLC.

Stability:

HBsAg Should be stored at 4°C, DO NOT FREEZE.