

# HBsAg adw Recombinant Protein



RPPB5523

---

## Product Information    Protein Information

**Product SKU:**

RPPB5523

**Host:**

Pichia Pastoris.

**Protein description:**

HbsAg adw produced Pichia Pastoris, having a molecular weight of approximately 24.0 kDa as shown on SDS-PAGE.

**Appearance:**

Sterile Filtered pale solution.

**Formulation:**

Sterile Filtered solution containing 20mM Phosphate Buffer, 154mM sodium chloride, pH 7.1.

**Purity:**

Greater than 95.0% as determined by SDS-PAGE.

**Stability:**

HBsAg Should be stored at 4°C.DO NOT FREEZE.