

CMV gB Recombinant Protein



RPPB5527

Product Information Protein Information

Product SKU:

RPPB5527

Host:

Escherichia Coli.

Protein description:

The E.Coli derived recombinant artificial mosaic protein contains the CMV gB immunodominant regions 11-67 amino acids, fused with a 26 kDa GST Tag, the protein weight is 6.5kDa, having a total weight of 32.5 kDa.

Formulation:

50mM Tris pH 7.2, 1mM EDTA and 50% glycerol.

Purity:

CMV gB protein is >95% pure as determined by 10% PAGE (coomassie staining).

Stability:

CMV gB protein although stable at 4°C for 1 week, should be stored below -18°C. Please prevent freeze thaw cycles.