

CMV Pp28 Recombinant Protein



RPPB5528

Product Information Protein Information

Product SKU:

RPPB5528

Host:

Escherichia Coli.

Protein description:

The E.Coli derived recombinant protein contains the CMV Pp28 (UL99) immunodominant regions, 130-160 amino acids.

Formulation:

25mM Tris-Hcl pH 7.2, 1mM EDTA and 50% glycerol.

Purity:

CMV Pp28 protein is >95% pure as determined by 10% PAGE (coomassie staining).

Stability:

CMV Pp28 protein although stable at 4°C for 1 week, should be stored below -18°C. Please prevent freeze thaw cycles.