

CMV Pp52 Recombinant Protein



RPPB5529

Product Information Protein Information

Product SKU:

RPPB5529

Host:

Escherichia Coli.

Protein description:

The E.coli derived 51 kDa recombinant protein contains the CMV Pp52 (UL44) immunodominant regions, 202-434 amino acids. Recombinant CMV-Pp52 is fused to a 26 kDa GST tag.

Formulation:

50mM Tris pH 7.2, 1mM EDTA and 50% glycerol.

Purity:

CMV Pp52 protein is >95% pure as determined by 10% PAGE (coomassie staining).

Stability:

CMV Pp52 protein although stable at 4°C for 1 week, should be stored below -18°C. Please prevent freeze thaw cycles.