

CMV Pp65 Recombinant Protein



RPPB5530

Product Information Protein Information

Product SKU:

RPPB5530

Host:

Escherichia Coli.

Protein description:

The E.Coli derived 52.2 kDa recombinant protein contains the CMV Pp65 (UL83) immunodominant regions, 297-510 amino acids. Recombinant CMV-Pp65 is fused to a 26 kDa GST tag.

Formulation:

CMV PP65 protein solution is supplied in PBS.

Purity:

CMV Pp65 protein is >90% pure as determined by SDS-PAGE.

Stability:

CMV Pp65 protein although stable at 4°C for 1 week, should be stored below -18°C. Please prevent freeze thaw cycles.