

CMV Pp65 Recombinant Protein



RPPB5531

Product Information Protein Information

Product SKU:

RPPB5531

Host:

Escherichia Coli.

Protein description:

The E.Coli derived 62.8 kDa recombinant protein contains the CMV Pp65 (UL83) immunodominant regions, having 561 amino acids.

Formulation:

25mM Tris-Hcl pH 8, 8M Urea, 5mM bMe.

Purity:

CMV Pp65 protein is >95% pure as determined by SDS-PAGE.

Stability:

CMV Pp65 protein although stable at 4°C for 1 week, should be stored below -18°C. Please prevent freeze thaw cycles.