

CMV Pp150 Recombinant Protein



RPPB5532

Product Information Protein Information

Product SKU:

RPPB5532

Host:

Escherichia Coli.

Protein description:

The E.Coli derived recombinant protein contains the CMV Pp150 (UL32) immunodominant regions, 1011-1048 amino acids.

Formulation:

50mM Tris pH 8.0, 1mM EDTA and 50% glycerol.

Purity:

CMV Pp150 protein is >95% pure as determined by 10% PAGE (coomassie staining).

Stability:

CMV Pp150 protein although stable at 4°C for 1 week, should be stored below -18°C. Please prevent freeze thaw cycles.