

CMV Mosaic Recombinant Protein



RPPB5533

Product Information Protein Information

Product SKU:

RPPB5533

Host:

Escherichia Coli.

Protein description:

The E.Coli derived Recombinant Cytomegalo Virus Mosaic contains the multiple epitopes from P150-P52-P38-P65-P28 having a molecular weight of 40kDa. The CMV Mosaic is fused to a 6xHis tag and purified by proprietary chromatographic technique.

Appearance:

Sterile filtered colorless solution.

Formulation:

CMV Mosaic protein solution (0.42mg/ml) contains Phosphate buffered saline, 2M urea and 100mM arginine.

Purity:

CMV Mosaic is >95% pure as determined by 12% PAGE (coomassie staining).

Stability:

CMV Mosaic protein although stable at 4°C for 1 week, should be stored below -18°C. Please prevent freeze thaw cycles.