

VZV gE Recombinant Protein



RPPB5534

Product Information Protein Information

Product SKU:

RPPB5534

Protein description:

The E.Coli derived recombinant protein contains the VZV gE immunodominant regions, amino acids 48-135 and fused to a GST-Tag at N-terminus.

Formulation:

60mM NaCl, 50mM Tris-HCl pH 8.0, 0.25% Sarkosil, 10mM glutathione and 50% glycerol.

Purity:

Varicella protein is >90% pure as determined by SDS-PAGE.

Stability:

Varicella Protein although stable at 4°C for 1 week, should be stored below -18°C. Please prevent freeze thaw cycles.