

VZV ORF9 Recombinant Protein



RPPB5535

Product Information Protein Information

Product SKU:

RPPB5535

Protein description:

The E.Coli derived recombinant protein contains the VZV ORF9 immunodominant regions, amino acids 6-28, 76-100.

Formulation:

25mM Tris-Hcl pH 8.0, 1mM EDTA and 50% glycerol.

Purity:

Varicella protein is >95% pure as determined by 10% PAGE (coomassie staining).

Stability:

Varicella protein although stable at 4°C for 1 week, should be stored below -18°C. Please prevent freeze thaw cycles.