

T.gondii MIC-3 Recombinant Protein



RPPB5537

Product Information Protein Information

Product SKU:

RPPB5537

Host:

Escherichia Coli.

Protein description:

The E.Coli derived recombinant protein contains the MIC3 immunodominant regions, amino acids 234-306. MIC3 protein is fused to a 26Kda Glutathione transferase tag.

Formulation:

10mM Tris-HCl pH 8, 1mM EDTA, 8M urea and 1mMDTT.

Purity:

Toxoplasma protein is >95% pure as determined by 10% PAGE (coomassie staining).

Stability:

Toxoplasma protein although stable at 4°C for 1 week, should be stored below -18°C. Please prevent freeze thaw cycles.