

T.gondii p24 Recombinant Protein



RPPB5538

Product Information Protein Information

Product SKU:

RPPB5538

Host:

Escherichia Coli.

Protein description:

The E.Coli derived recombinant protein contains the p24 (GRA1) immunodominant regions, amino acids 57-149.

Formulation:

(1mg/ml) 50mM Tris pH 8, 1.5M urea and 50% glycerol.

Purity:

Toxoplasma protein is >95% pure as determined by 10% PAGE (coomassie staining).

Stability:

Toxoplasma Protein although stable at 4°C for 1 week, should be stored below -18°C. Please prevent freeze thaw cycles.