

Toxoplasma P22 Recombinant Protein



RPPB5539

Product Information Protein Information

Product SKU:

RPPB5539

Host:

Escherichia Coli.

Protein description:

Recombinant Toxoplasma Gondii P22 (SAG2) produced in E.coli is a full length Surface antigen 2 (SAG2) and fused with a 6xHis tag at C- terminus.

Appearance:

Sterile Filtered solution.

Formulation:

Sterile Filtered solution containing Phosphate buffered saline.

Purity:

Toxoplasma protein is >95% pure as determined by 12% PAGE (coomassie staining).

Stability:

Toxoplasma P22 although stable at 4°C for 1 week, should be stored below -18°C. Please prevent freeze thaw cycles.