

# T.gondii p29 Recombinant Protein



RPPB5540

---

## Product Information    Protein Information

**Product SKU:**

RPPB5540

**Host:**

Escherichia Coli.

**Protein description:**

The E.coli derived recombinant protein contains the p29 (GRA7) immunodominant regions, amino acids 24-100 and fused to a 26kDa GST tag.

**Formulation:**

50mM Tris-HCl, 1.5M urea and 50% glycerol.

**Purity:**

Toxoplasma protein is >95% pure as determined by 10% PAGE (coomassie staining).

**Stability:**

Toxoplasma protein although stable at 4°C for 1 week, should be stored below -18°C. Please prevent freeze thaw cycles.