

# T.gondii p30 Recombinant Protein



RPPB5541

---

## Product Information    Protein Information

**Product SKU:**

RPPB5541

**Host:**

Escherichia Coli.

**Protein description:**

The E.Coli derived recombinant protein contains the p30 (SAG1) immunodominant regions, amino acids 45-198 and fused to a 26kDa GST tag, having a total Mw of 42.3 kda.

**Formulation:**

50mM Tris-HCl pH 8, 1.5M urea and 50% glycerol.

**Purity:**

Toxoplasma protein is >95% pure as determined by 10% PAGE (coomassie staining).

**Stability:**

Toxoplasma protein although stable at 4°C for 1 week, should be stored below -18°C. Please prevent freeze thaw cycles.