

Toxoplasma P32 Recombinant Protein



RPPB5542

Product Information Protein Information

Product SKU:

RPPB5542

Host:

Escherichia Coli.

Protein description:

Recombinant Toxoplasma Gondii P32 (GRA6) produced in E.coli containing 180 amino acids of the Dense Granule Antigen 6 (GRA6) and fused with a 6xHis tag at C- terminus. This antigen forms dimers on SDS-PAGE under denatured conditions, compared to other toxoplasma monomer antigens, it has better immunoreactions.

Appearance:

Sterile Filtered solution.

Formulation:

Sterile Filtered solution containing containing PBS and 25mM arginine.

Purity:

Toxoplasma protein is >95% pure as determined by 10% PAGE (coomassie staining).

Stability:

Toxoplasma P32 although stable at 4°C for 1 week, should be stored below -18°C. Please prevent freeze thaw cycles.