

T.gondii p40 Recombinant Protein



RPPB5543

Product Information Protein Information

Product SKU:

RPPB5543

Host:

Escherichia Coli.

Protein description:

The E.Coli derived recombinant Toxoplasma Gondii p40 having a Mw of 35kDa and fused to a 6 His tag at C-terminus can be used to test the specific IgG and IgM antibody for the diagnosis of Toxoplasma gondii infection.

Appearance:

Sterile Filtered clear solution.

Formulation:

1x PBS and 25mM arginine.

Purity:

Toxoplasma protein is >95% pure as determined by 10% PAGE (coomassie staining).

Stability:

T.gondii p40 although stable at 4°C for 1 week, should be stored below -18°C. Please prevent freeze thaw cycles.